Retreat Platinum Water Jet 60" x 60" x 23" **.** O

SKU: RETREAT-PWJ

TECHNICAL SPECIFICATIONS		BA	THERS: 2	
Drain	Side drain	INSTALLATIO		
Pump Location(s): Select 2	Std. = Standard / Alt LF/RF Optional	•	DECK (D)	
Bath Weight (empty)	180 lbs		UNDER MO	
Gallons to Operate	50-75		ALCOVE (A	
ELECTRICAL SPECIFICATIONS			FREESTAN	
Motor/Pump	(2) 115v 20AMP GFCI Circuits		CLAWFOOT	
Inline Stainless Steel Heater	(1) 115v 20AMP GFCI Circuit		PEDESTAL	

STANDARD FEATURES & ACCESSORIES for Retreat Platinum Water Jet

JETS

15 Water Jets Back Massage System(s) (7 Jets) Neck Massage Pillow (2 Jets) Sole Soother Jet(s) 🛞 Accumassage Back, Neck⊕ Ultra Flow Jets -

PUMPS

One - 3.25 & One - 2.25 Hydropower 3-speed Pump Quiet Pad Anti-Vibration Mounted Pump

HEATERS

Stainless Steel Electric Heater

CONTROLS

Air Switch On/Off Control(s) For Chromatherapy Lighting Air Switch On/Off with digital Temperature Readout Control 3 Speed Massage: High, Medium or Low

PURITY

HydroClean Water Purity System HydroCleanse

INSTALLATION DECK (D) UNDER MOUNT (U) ALCOVE (A) FREESTANDING (F) CLAWFOOT (L) PEDESTAL (P)





STANDARD ACCESSORIES


Grab Bar(s) - 1 on Side Drain baths White, Biscuit, Bone, Brushed Nickel, Polished Chrome Chromatherapy Lighting Deluxe Cable Drain - White, Biscuit, Bone, Brushed Nickel, Oil Rubbed Bronze, Polished Chrome Deluxe Cushion Pillow on two person baths Neck Massage Pillow EZ Level Base FDA Recipe Non-Toxic Food & Beverage Grade Plumbing (primarily) 30 Point Quality Assurance Factory Test

REQUIREMENTS

Pump Locations: Select 2 Electrical: three 115v 20AMP GFCI Circuits Drain and Grab Bars: Select Color / Finish Tub should not be supported by the lip in any way. Level tub & screw or glue EZ Level adjustable legs to sub floor. CUT OUT DIMENSIONS: **OVALS:** A template will be provided

RECTANGLES: Subtract 1" to 1.5" from length & width for hole cut out size

EZ CLEAN CAST ACRYLIC



NOTE: Graphics may not be to scale. Specifications and features are subject to change without notice and are ± .25". All final measurements should be taken from tub on site. Access 18"x 18" or larger is required for each pump. Some pump locations and plumbing may extend beyond the perimeter of the bath.