#### Zen 6632 SD Platinum Combination 66" x 32" x 23" ± 000

SKU: ZEN6632SD-PCB

TECHNICAL SPECIFICATION	IS	
Drain	Side drain	
Pump Location(s): Select 3	Std. = Standard / Alt LF/RF Optional	▎▕
Bath Weight (empty)	155 lbs	] [
Gallons to Operate	44-69	] [
ELECTRICAL SPECIFICATIONS		
Motor/Pump	(3) 115v 20AMP GFCI Circuits	] [
Inline Stainless Steel Heater	(1) 115v 20AMP GFCI Circuit	] [

### **STANDARD FEATURES & ACCESSORIES for PLATINUM COMBINATION**

#### JETS

38 Air & Water Massage Jets Back Massage System(s) (7 Jets) Neck Massage Pillow (2 Jets) ( $\oplus \oplus$ ) Sole Soother Jet(s) Accumassage Back & Neck Jets  $\oplus$ Ultra Flow Jets -Classic Air Jets O EZ Clean Air Jets (optional) X

#### PUMPS

Two-3.25 Hydropower 3-speed Pumps (Side Drain Tubs) 1.25 Air Power Variable Speed, Pulse & Wave heated Air Pump Quiet Pad Anti-Vibration Mounted Pump

#### **HEATERS**

Stainless Steel Heater Warm Back Rest Tubing Zone

#### CONTROLS

Air Switch On/Off for Chromatherapy Lighting Air Switch On/Off with digital Temperature Readout Hydro Pumps Electronic Variable Speed Air Massage Control with Automatic Purge & Dry Cycle Control Air Pump 3 Speed Massage: High, Medium or Low

#### PURITY

HydroClean Water Purity System HydroCleanse

## BATHERS: 1 INSTALLATION DECK (D) UNDER MOUNT (U) ALCOVE (A)

FREESTANDING (F) CLAWFOOT (L) PEDESTAL (P)

Jet & Control placement may vary

# PLATINUM COMBINATION



#### STANDARD ACCESSORIES

Grab Bar(s) - 1 on Side Drain Baths White, Biscuit, Bone, Brushed Nickel, Polished Chrome Chromatherapy Lighting Deluxe Cable Drain - White, Biscuit, Bone, Brushed Nickel, Oil Rubbed Bronze, Polished Chrome Deluxe Cushion Pillow on 2 person baths Neck Massage Pillow EZ Level Base FDA Recipe Non-Toxic Food & Beverage Grade Plumbing (primarily) 30 Point Quality Assurance Factory Test


#### REQUIREMENTS

Pump Locations: Select 3 Electrical: Four 115v 20AMP GFCI Circuits Cable Drain and Grab Bars: Select Color / Finish Tub should not be supported by the lip in any way. Level tub & screw or glue adjustable EZ Level legs to sub floor. **CUT OUT DIMENSIONS:** OVALS: A template will be provided

RECTANGLES: Subtract 1" to 1.5" from length & width for hole cut out size



- A K HTT® Leg Therapy Air Massage System includes 10 Air Jets on the floor. Optional on Combo Plat and Combo Gold Series
- **B** X 6 Jet HTT® Back Massage System: N/A on Combo Series-Standard on Plat series



Some standard & most alternate pump locations may stickout past perimeter of the bath 1"-3". NOTE: Graphics may not be to scale. Specifications and features are subject to change without notice and are ± .25". All final measurements should be taken from tub on site. Access 18"x 18" or larger is required for each pump. Consult your contractor for installation requirements.